

# AWI-Loc 1.5



## **206075**

*Attachment Matritze grün T4®-switch / plus (VE = 4)*

*Divergenzausgleich: 20 – 40°*

*Abzugskrafttrange: 23,5 – 27,5 N*

*Attachment Matrix green T4®-switch / plus (SP = 4)*

*divergence compensation: 20 – 40°*

*reinforced retention: 23.5 – 27.5 N*



## **206077**

*Attachment Matritze orange T4®-switch / plus (VE = 4)*

*Divergenzausgleich: 20 – 40°*

*Abzugskrafttrange: 13 – 16,5 N*

*Attachment Matrix orange T4®-switch / plus (SP = 4)*

*divergence compensation: 20 – 40°*

*reinforced retention: 13 – 16.5 N*



## **206076**

*Attachment Matritze rot T4®-switch / plus (VE = 4)*

*Divergenzausgleich: 20 – 40°*

*Abzugskrafttrange: 5 – 8 N*

*Attachment Matrix red T4®-switch / plus (SP = 4)*

*divergence compensation: 20 – 40°*

*reinforced retention: 5 – 8 N*



## **206078**

*Titan-Hülse verbaut mit Matritze schwarz*

*T4®-switch / plus*

*Titanium-Sleeve incl. Matrix black*

*T4®-switch / plus*

# AWI-Loc 1.5



## **206083**

*Matritzenset-2 T4®-switch / plus inkl.*

*Block-out-Ring*

*Divergenzausgleich: 20 – 40°*

*Matrix Set-2 T4®-switch / plus incl. Block-out-Ring*

*divergence compensation: 20 – 40°*



## **206081**

*Laborimplantat für ClicLoc T4®-switch / plus*

*Implant analog ClicLoc T4®-switch / plus*



## **206080**

*Abdruckpfosten verbaut mit Matritze schwarz*

*T4®-switch / plus*

*Impression post incl. Matrix black*

*T4®-switch / plus*



## **206084**

*Universal-Einbringinstrument ClicLoc*

*T4®-switch / plus*

*Universaltool ClicLoc T4®-switch / plus*